

\\USER\FORSKNINGSPROJEKT\NEURO\BoF 127 MM\FWF_MENI_MAX_STE

TA: 4:12 PM: FIX Voxel size: 2.3×2.3×2.3 mmPAT: 2 Rel. SNR: 1.00 : epse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	40
Dist. factor	0 %
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
TR	6000 ms
TE	80,0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HE1-4

Contrast - Common

TR	6000 ms
TE	80,0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
Base resolution	100
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	EPI/separate

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off
Dynamic Field Corr.	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	40
Dist. factor	0 %
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
TR	6000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.2 A18.7 H5.4
R	2,2 mm
A	18,7 mm
H	5,4 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T

System - Miscellaneous

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Rotation	0,00 deg
A >> P	230 mm
R >> L	230 mm
F >> H	92 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123,259684 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

Physio - Signal1

1st Signal/Mode	None
TR	6000 ms
Concatenations	1

Physio - PACE

Resp. control	Off
Concatenations	1

Diff - Neuro

Diffusion mode	Free
Diff. directions	39
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	2000 s/mm ²
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	40

Diff - Body

Diffusion mode	Free
Diff. directions	39
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	2000 s/mm ²
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm ²
Noise level	40

Diff - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0,93 ms
Bandwidth	1190 Hz/Px

Sequence - Part 2

EPI factor	100
RF pulse type	Normal
Gradient mode	Fast
Excitation	Standard

Sequence - pTX Pulses**Sequence - Special**

WFBank	Bank06
WFSelect	Cust01
RotationMode	XYZ
NormalizeMode	To max
PostWFMode	RevAmp
TimingMode	Manual
PauseMode	Min
BalanceGradMode	On
ICEHeader	Standard
MaxBVal	2162 s/mm ²
PreDur	29500 µs
PostDur	21240 µs
PauseDur	8020 µs

\\USER\FORSKNINGSPROJEKT\NEURO\BoF 127 MM\FWF_MENI_MAX_LTE

TA: 4:12 PM: FIX Voxel size: 2.3×2.3×2.3 mmPAT: 2 Rel. SNR: 1.00 : epse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	40
Dist. factor	0 %
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
TR	6000 ms
TE	80,0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HE1-4

Contrast - Common

TR	6000 ms
TE	80,0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
Base resolution	100
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
----------	--------

Resolution - iPAT

Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	EPI/separate

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off
Dynamic Field Corr.	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	40
Dist. factor	0 %
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
TR	6000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.2 A18.7 H5.4
R	2,2 mm
A	18,7 mm
H	5,4 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T

System - Miscellaneous

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Rotation	0,00 deg
A >> P	230 mm
R >> L	230 mm
F >> H	92 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123,259684 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

Physio - Signal1

1st Signal/Mode	None
TR	6000 ms
Concatenations	1

Physio - PACE

Resp. control	Off
Concatenations	1

Diff - Neuro

Diffusion mode	Free
Diff. directions	39
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	2000 s/mm ²
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	40

Diff - Body

Diffusion mode	Free
Diff. directions	39
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	2000 s/mm ²
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm ²
Noise level	40

Diff - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0,93 ms
Bandwidth	1190 Hz/Px

Sequence - Part 2

EPI factor	100
RF pulse type	Normal
Gradient mode	Fast
Excitation	Standard

Sequence - pTX Pulses**Sequence - Special**

WFBank	Bank06
WFSelect	Cust02
RotationMode	XCh-1D
NormalizeMode	To max
PostWFMode	None
TimingMode	Equal
PauseMode	Min
BalanceGradMode	On
ICEHeader	Standard
MaxBVal	5665 s/mm ²
PreDur	23620 μ s
PostDur	23620 μ s
PauseDur	8020 μ s

\\USER\FORSKNINGSPROJEKT\NEURO\BoF 127 MM\FWF_MENI_MAX_STE_PA

TA: 0:36 PM: FIX Voxel size: 2.3×2.3×2.3 mmPAT: 2 Rel. SNR: 1.00 : epse

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	40
Dist. factor	0 %
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
TR	6000 ms
TE	80,0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HE1-4

Contrast - Common

TR	6000 ms
TE	80,0 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms

Resolution - Common

FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
Base resolution	100
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2

Resolution - iPAT

Ref. lines PE	24
Reference scan mode	EPI/separate

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off
Dynamic Field Corr.	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	40
Dist. factor	0 %
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	230 mm
FoV phase	100,0 %
Slice thickness	2,3 mm
TR	6000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.2 A18.7 H5.4
R	2,2 mm
A	18,7 mm
H	5,4 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

Geometry - Saturation

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

System - Miscellaneous

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R2.2 A18.7 H5.4 mm
Orientation	Transversal
Rotation	0,00 deg
A >> P	230 mm
R >> L	230 mm
F >> H	92 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123,259684 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

Physio - Signal1

1st Signal/Mode	None
TR	6000 ms
Concatenations	1

Physio - PACE

Resp. control	Off
Concatenations	1

Diff - Neuro

Diffusion mode	Orthogonal
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	100 s/mm ²
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	40

Diff - Body

Diffusion mode	Orthogonal
----------------	------------

Diff - Body

Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	100 s/mm ²
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm ²
Noise level	40

Diff - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0,93 ms
Bandwidth	1190 Hz/Px

Sequence - Part 2

EPI factor	100
RF pulse type	Normal
Gradient mode	Fast
Excitation	Standard

Sequence - pTX Pulses**Sequence - Special**

WFBank	Bank06
WFSelect	Cust01
RotationMode	XYZ
NormalizeMode	To max
PostWFMode	RevAmp
TimingMode	Manual
PauseMode	Min
BalanceGradMode	On
ICEHeader	Standard
MaxBVal	2162 s/mm ²
PreDur	29500 μ s
PostDur	21240 μ s
PauseDur	8020 μ s